


# HPV ONCOTECT™ E6, E7 Kit

## Quantify and Monitor HPV Gene Expression



### ULTRA SENSITIVITY

- High throughput Cervical Cancer Screening using Intracellular HPV E6 & E7 mRNA Quantification and Flow Cytometry
- Detects HPV mRNA at 5-10 copies per cell.

### ULTRA SPECIFICITY

- HPV E6 & E7 are consistently expressed in cervical cancer cell lines, HPV transformed cells and in HPV-associated pre-cancer and cancer.

### VERSATILE


- Used with PAP Smears, fine-needle aspirates, punch biopsies, cell lines, cell suspensions and tissue samples.
- Cellular morphology preserved for cytopathological evaluation.

### SIMPLE-FAST RESULTS

- Assay Time - 2 Hours compared to 2 days for typical PCR or hybrid capture assays determining serotype only.
- Can be run on any flow cytometry platform.

Kits include all reagents and FL\*-labeled HPV probes

- 100-Test Kit, Model 88002



Fluorescent micrograph of cervical sample demonstrating HPV E6 & E7 mRNA processed with HPV OncoTect.

For More Information Contact: [www.invirion.com](http://www.invirion.com)  
Invirion Diagnostics, LLC  
1900 Spring Road, Suite 200 Oak Brook Illinois 60523  
Phone: 630-954-2875

  
Diagnostics, LLC